

Atlantis® (IO) FLO-S compatibility chart

The scan body Atlantis (IO) FLO-S is used when scanning for Atlantis suprastructures. It can be used on implant/abutment level for all positions in the mouth of the patient or on the master model during scanning. Refer to instructions for use for further information, 2663-IFU Atlantis IO FLO-S.

This chart describes what implants/abutments the Atlantis (IO) FLO-S is compatible with and the content of Dentsply Sirona Implants Atlantis suprastructures (IO) FLO-S libraries. The libraries are used for digital design guidance for Atlantis suprastructures and for building and printing models. The libraries can be used for 3Shape and Exocad design software.

Note: The FLO-S scan body without the IO marking can only be used for model scanning and is not approved for intraoral scanning. The IO FLO-S scan body will be used for both intraoral and model scanning.

Each Atlantis (IO) FLO-S reference number is sold as a kit containing a tray with 6 (IO) FLO-S for storing the non-sterile (IO) FLO-S and a Tag (to be placed in the Atlantis FLO-S case (REF) 36701) for 1 specific connection.



Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Dentsply Sirona	IMPLANT LEVEL							
	EV Connection: PrimeTaper EV Implant / OmniTaper EV Implant / Astra Tech Implant EV							
	EV Connection (XS) Green	A04K	68020018	YES	YES	PMA-ATE30-1		
	EV Connection (S) Purple	A04L	68020019	YES	YES	PMA-ATE36-1		
	EV Connection (M) Yellow	A04M	68020020	YES	YES	PMA-ATE42-1		
	EV Connection (L) Blue	A04N	68020021	YES	YES	PMA-ATE48-1		
	EV Connection (XL) Brown*	A04O	68020022	YES	YES	PMA-ATE54-1		
	Astra Tech Implant EV Profile							
	4.2 Yellow	A04M	68020020	YES	YES	PMA-ATE42P-1		
	4.8 Blue	A04N	68020021	YES	YES	PMA-ATE48P-1		
	Astra Tech Implant System TX							
	OsseoSpeed TX 3.0 – Yellow	A04C	68020023	YES	YES	PMA-ATO30-1		
	OsseoSpeed TX 3.5, 4.0 – Aqua	A04D	68020024	YES	YES	PMA-ATO35-1		
	OsseoSpeed TX 4.5, 5.0 – Lilac	A04E	68020025	YES	YES	PMA-ATO45-1		
	OsseoSpeed TX Profile 4.5/5.0	A04E	36741-6		YES			
	Xive / Frialit							
	Xive S 3.0	F01A	68020026	YES	YES	PMA-DXI30-1		
	Xive S 3.4 / Frialit 3.4	F01B	68020027	YES	YES	PMA-DXI34-1		
	Xive S 3.8 / Frialit 3.8	F01C	68020028	YES	YES	PMA-DXI38-1		
	Xive S 4.5 / Frialit 4.5	F01D	68020029	YES	YES	PMA-DXI45-1		
	Xive S 5.5 / Frialit 5.5	F01E	68020030	YES	YES	PMA-DXI55-1		

* **Note:** Not applicable for PrimeTaper EV Implant.

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

Note: Atlantis (IO) FLO-S assortment can vary from country to country. Please check availability with your local Customer Service Representative.

Please turn over ➡

Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information					
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging		
Dentsply Sirona	ABUTMENT LEVEL									
	EV Connection: DS PrimeTaper Implant / DS OmniTaper Implant / Astra Tech Implant EV									
	MultiBase Abutment EV	A04S	68020034	YES	YES	PMA-MUA45-1	HBN-MUA45-1	HBN-MUA45H-1		
	Uni Abutment EV	A04F	68020031	YES	YES	PMA-AUE33-1	HBN-AUE33-1	HBN-AUE33H-1		
	Astra Tech Implant System TX									
	20° UniAbutment	A04A	68020032	YES	YES	PMA-ATU20-1	HBN-ATU20-1	HBN-ATU20H-1		
	45° UniAbutment	A04B	68020033	YES	YES	PMA-ATU45-1	HBN-ATU45-1			
	Xive									
	Xive MP/TG 3.4,3.8	F01F	36745-6		YES					
	Xive MP/TG 4.5	F01G	36746-6		YES					
	Xive MP 5.5	F01H	36747-6		YES					
	Ankylos									
	Ankylos Balance Base Narrow D4.2	A03A	68020035	YES	YES	PMA-DBA40-1	HBN-DBA40-1	HBN-DBA40H-1		
	Ankylos Balance Base D5.5	A03B	36730-6		YES					
MIS	IMPLANT LEVEL									
	Internal Hex Narrow SEVEN - M4	M04D	68020043	YES	YES					
	Internal Hex Standard SEVEN - M4	M04A	68020044	YES	YES					
	Internal Hex Wide SEVEN - M4	M04B	68020045	YES	YES					
	Conical Connection Narrow C1	M04E	68020046	YES	YES					
	Conical Connection Standard C1 - V3	M04F	68020047	YES	YES					
	Conical Connection Narrow V3	M04E	68020046	YES	YES					
	Conical Connection Wide C1	M04G	68020048	YES	YES					
	ABUTMENT LEVEL									
	Multi-Unit Abutment	M04H	68020049	YES	YES					
	CONNECT abutment D4	M04I	68021007							
	CONNECT abutment D5.7	M04J	68021008							
	BioHorizons	IMPLANT LEVEL								
		Internal/ Tapered 3.5	B03A	36710-6						
Internal/ Tapered 4.5		B03B	36711-6							
Internal/ Tapered 5.7		B03C	36712-6							
ABUTMENT LEVEL										
Multi-unit abutment		A04S	68020034							

Please turn over ➡

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

Note: Atlantis (IO) FLO-S assortment can vary from country to country.
Please check availability with your local Customer Service Representative.



Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Biomet 3i	IMPLANT LEVEL							
	Certain 3.25, 4/3	301F	36716-6					
	Certain 4.0, 5/4	301H	36717-6					
	Certain 5.0	301J	36718-6					
	Certain 6.0	301L	36719-6					
	Certain Prevail 3/4/3, 4/3	301F	36716-6					
	Certain Prevail 4/5/4, 5/4	301H	36717-6					
	Certain Prevail 5/6/5, 6/5	301J	36718-6					
	Certain XP4/5	301J	36718-6					
	Certain XP 5/6	301L	36719-6					
	MicroMiniplant 3.25 / Miniplant 3.25	301G	36721-6					
	Osseotite 3.75 / 4.0	N01F	68020009					
	Osseotite 5.0, XP 4/5	301K	36723-6					
	Osseotite 6.0, XP 5/6	301M	36724-6					
	ABUTMENT LEVEL							
Low Profile Abutment	301E	36720-6						
Camlog	IMPLANT LEVEL							
	Screw-line Implant 3.3	C01E	36725-6					
	Screw-line/Root-line Implant 3.8	C01G	36726-6					
	Screw-line/Root-line Implant 4.3	C01H	36727-6					
	Screw-line/Root-line Implant 5.0	C01I	36728-6					
	Screw-line/Root-line Implant 6.0	C01K	36729-6					
Keystone Dental	IMPLANT LEVEL							
	PrimaConnex SD 3.3,3.5	K01A	36755-6					
	PrimaConnex RD 4.0,4.1	K01B	36754-6					
	PrimaConnex WD 5.0	K01C	36756-6					
	Genesis 3.8	K01A	36755-6					
	Genesis 4.5	K01B	36754-6					
Genesis 5.5, 6.5	K01C	36756-6						
Megagen	IMPLANT LEVEL							
	AnyOne 3.5	M03A	68021002					
	AnyOne 4.0, 4.5, 5.0, 5.5, 6.0, 6.5, 7.0, 7.5, 8.0	M03A	68021002					
	ABUTMENT LEVEL							
	Multi-unit Abutment (N_Type)	A04S	68020034					
	Multi-unit Abutment (S_Type)	M03G	68021004					
	Octa Abutment 3.8	M03H	68021005					
	Octa Abutment 4.8	S05C	36772-6					
	Octa Abutment 5.8	M03I	68021006					

Please turn over 

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

Note: Atlantis (IO) FLO-S assortment can vary from country to country.
Please check availability with your local Customer Service Representative.



Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Neodent	IMPLANT LEVEL							
	GM Implant 3.5, 3.75, 4.0, 4.3	N03A	68021009					
	GM Implant 5.0, 6.0, 7.0	N03A	68021009					
	ABUTMENT LEVEL							
	GM Abutment	N03C	68021010					
	GM Micro Abutment	N03B	68021011					
	GM Mini Conical Abutment	N01H	36759-6					
Nobel Biocare	IMPLANT LEVEL							
	NobelActive NP 3.5	N01A	68020050	YES	YES	PMA-NBA35-1		
	NobelActive RP 4.3, 5.0	N01B	68020051	YES	YES	PMA-NBA43-1		
	NobelActive WP 5.5	N01O	68020052	YES	YES	PMA-NBA60-1		
	NobelReplace NP 3.5	N01J	68020053	YES	YES	PMA-NBR35-1		
	NobelReplace RP 4.3	N01K	68020054	YES	YES	PMA-NBR43-1		
	NobelReplace WP 5.0	N01L	68020055	YES	YES	PMA-NBR50-1		
	NobelReplace 6.0	N01M	68020056	YES	YES	PMA-NBR60-1		
	Branemark System NP 3.3	N01E	68020057	YES	YES	PMA-BRA35-1		
	Branemark System RP 3.75, 4.0	N01F	68020009	YES	YES	PMA-BRA41-1		
	ABUTMENT LEVEL							
Multi-Unit Abutment RP	N01H	68020058	YES	YES	PMA-MUA45-1			
Straumann	IMPLANT LEVEL							
	Bone Level 3.3 NC	S05G	68020059	YES	YES	PMA-SBO33-1		
	Bone Level 4.1, 4.8 RC	S05H	68020060	YES	YES	PMA-SBO41-1		
	Standard/Standard Plus 4.8 RN (Regular Neck)	S05D	68020061	YES	YES	PMA-SSY48-1		
	Standard/Standard Plus 4.8 WN (Wide Neck)	S05F	68020062	YES	YES	PMA-SSY65-1		
	BLX Implant RB 3.75, 4.5	S05S	68021015					
	BLX Implant WB 5.5, 6.5	S05S	68021015					
	ABUTMENT LEVEL							
	Monotype 8° Conus	S05A	36771-6		YES			
	RN Abutment Level	S05C	36772-6		YES			
	WN Abutment Level	S05E	36777-6		YES			
	Bone Level Multi-Base Abutment D3.5	S05I	36768-6		YES			
	Bone Level Multi-Base Abutment D4.5	S05J	36769-6		YES			
	Bone Level Multi-Base Angled Abutment	S05K	36770-6		YES			
	Screw-Retained Abutment 3.5	S05O	68020069	YES	YES	PMA-SSA35-1		
Screw-Retained Abutment 4.6	S05P	68020070	YES	YES	PMA-SSA46-1			

Please turn over ➡

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

Note: Atlantis (IO) FLO-S assortment can vary from country to country.
Please check availability with your local Customer Service Representative.



Manufacturer	Implant system	FLO-S id	FLO-S Movex no (Kit)	Intraoral scanning	Library information			
					Available in library	Elos Accurate Analog for Printed Models	Elos Accurate Hybrid Base Non-Engaging	Elos Accurate Hybrid Base H Non-Engaging
Zimmer Dental	IMPLANT LEVEL							
	Tapered Screw-Vent™ 3.5 / Screw Vent™ 3.3, 3.7	B03A	36710-6					
	Tapered Screw-Vent™ 4.5 / Screw Vent™ 4.5	B03B	36711-6					
	Tapered Screw-Vent™ 5.7	B03C	36712-6					
	ABUTMENT LEVEL							
	Tapered Abutment	Z01E	36779-6					
	Spectra Cone Abutment	Z01F	36778-6					

All trademarks and company names are the property of their respective owners.

Elos Accurate® is a registered trademark.

Note: Atlantis (IO) FLO-S assortment can vary from country to country.
Please check availability with your local Customer Service Representative.